# BIGINELLI DIHYDROPYRIMIDINE SYNTHESIS

Reflux and Thin Layer Chromatography (TLC)

## **EXPERIMENTAL OBJECTIVES**

To synthesize an organic compound and examine its purity using thin layer chromatography (TLC)

## **LEARNING OBJECTIVES**

- To learn how to carry a reaction from start to finish, including reaction setup, monitoring, workup, purification, and product verification
- To learn how to setup and why to use a reflux reaction
- To learn how to use TLC to monitor a reaction and determine sample purity

#### PRE-LAB EXERCISES

Please read:

See links on Blackboard for **Reflux**, and **Thin-layer Chromatography**Calculate amounts of starting materials by converting from millimoles (mmol) to grams.

#### **BACKGROUND**

The Biginelli pyrimidone synthesis is a reaction named after the scientist who first published its use in 1891, Pietro Biginelli.<sup>1</sup>

The multicomponent reaction involves the combining of a  $\beta$ -dicarbonyl compound (also known as a 1,3-dicarbonyl), an aldehyde, and a urea to form a dihydropyrimidine heterocycle. The reaction mechanism is quite complex, involving an aldol condensation (Carey Chpt. 20), a 1,4-addition (Carey Chpt. 20), and subsequent ring closure with elimination (Carey Chpt 5). The dihydropyrimidine heterocycle product is a common pharmacophore and is found in molecules with antiviral, antitumor, antibacterial, anti-inflammatory, and antihypertensive activity.<sup>2</sup> For

<sup>&</sup>lt;sup>1</sup> Biginelli, P. Aldehyde-Urea Derivatives of Aceto- and Oxaloacetic Acids. Gazz. Chim. Ital. 1893, 23, 360-413

<sup>&</sup>lt;sup>2</sup> Kappe, C. O. Recent Advances in the Biginelli Dihydropyrimidine Synthesis. Accounts of Chemical Research. 2000, 33, 879-888.

our purposes though, it is a very efficient reaction, which can be monitored by TLC, and it produces an easily isolatable product. Specifically, we will react 4-methoxybenzaldehyde, urea, and ethylacetoacetate (ethyl 3-oxobutanoate) under acid catalysis to form ethyl 4-(4-methoxyphenyl)-6-methyl-2-oxo-1,2,3,4-tetrahydropyrimidine-5-carboxylate.

## PROCEDURE - WEEK 1

- 1. Place 10 mmol vanillin, 15 mmol ethyl acetoacetate, 10 mmol urea, and 5 mL of 95% ethanol into a 25 mL round bottom flask equipped with a reflux condenser
- 2. Add 0.2 mL of concentrated HCl to the reaction solution
- 3. Heat the reaction mixture at reflux for 1.5 hours
- 4. Cool the flask to 0°C, a precipitate should form
- 5. Collect the solid product by vacuum filtration, wash with cold 95% ethanol (~5 mL)
- 6. Once dry, be sure to record the mass of your solid product and save for the next lab period
- 7. Discard all waste into the appropriate waste container

## PROCEDURE - WEEK 2

- 1. Analyze your product using Thin Layer Chromatography (TLC).
  - a. Set up a co-spot TLC plate using vanillin and your crude reaction product
  - b. Run the TLC in a closed chamber using your chosen mobile phase
  - c. When the solvent line nears the top of the plate, remove the plate from the chamber, mark the solvent line, and then allow it to dry
  - d. Visualize your spots using a UV-light (be sure to outline the spots in pencil on the TLC plate)
  - e. Determine  $R_{\rm f}\, values$  for each spot, and the purity of your crude product
  - f. Copy the TLC plate into your lab notebook and discard the TLC plate into the broken glass container
- 2. After TLC analysis, make a determination of how pure your product is...
- 3. Also verify the purity and identity of your compound by melting point analysis (your need to report a range expected melting point ~220°C).
- 4. Discard your product in the appropriate waste container when you are certain you have acquired all necessary data.